**erbe**  
power your performance.

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**HybridKnife I-Typ**  
Ø 2,3 mm; Länge 1,9 m  
Nr. 20150-061

**HybridKnife O-Typ**  
Ø 2,3 mm; Länge 1,9 m  
Nr. 20150-062

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## Mukosa-Elevation (Flexible Sonde + HybridKnife):

Ösophagus/Magen	ERBEJET 2, Effekt 30–50
Rechtes Kolon	ERBEJET 2, Effekt 10–15
Rektum/linkes Kolon	ERBEJET 2, Effekt 20–30

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## Referenzen

*Repici, A et al.: A prospective, single center study of endoscopic submucosal dissection of rectal LST lesions larger than 3 cm by using an innovative concept of injecting and Cutting: The water-Jet Hybrid-Knife (ESD-H). Gastrointest Endosc Vol 73, Issue 4, Supplement, Page AB156, April 2011*

*Fernandez-Esparrach, G. et al.: A novel device for endoscopic submucosal dissection that combines water-jet submucosal hydrodissection and elevation with electrocautery: initial experience in a porcine model. Gastrointest Endosc Vol 71, No. 3: 2010*

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*Neuhaus, H.: Endoscopic submucosal dissection in the upper gastrointestinal tract: Present and Future View of Europe. Digestive Endoscopy (2009) 21 (Suppl. 1), S4–S6*

*Neuhaus, H.: Endoscopic mucosal resection and endoscopic submucosal dissection in the West – too many concerns and caveats? Endoscopy 2010; 42: 859–861*

*Schumacher, B. et al.: Endoscopic submucosal dissection of early gastric neoplasia with a water jet–assisted knife: a Western, single-center experience. Gastrointest Endosc, 2012.*

*Der direkte Link zur EMR  
[www.medical-video.com](http://www.medical-video.com)*

*Der direkte Link zur ESD  
[www.medical-video.com](http://www.medical-video.com)*

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Erbe Elektromedizin GmbH