



Instruments and accessories

for hydrosurgery



Legend


Important notes

The authors and Erbe Elektromedizin GmbH have taken the greatest possible care while compiling the catalogue with electrosurgical accessories. Nevertheless, it is impossible to exclude the possibility that errors may be contained herein.

The information and recommendations given here may not be construed as constituting any basis for any claims against Erbe Elektromedizin GmbH. Should legal regulations stipulate liability, such liability shall be limited to intentional misconduct or gross negligence.

The data on the use of instruments is based on clinical experience. Medical science is subject to a process of permanent change based on research and clinical experience. Therefore the products are subject to alterations in design and technique. Because of production tolerances dimensions in the catalogue may marginally differ.

Further product information



Please read carefully the corresponding notes on use before using any instrument.

Hydrosurgery

Instruments and accessories

- 04 Hybrid instruments
- 06 Applicators and probes
- 08 Accessories for ERBEJET 2
- 09 Accessories for ESM 2
- 11 Carts and accessories

Instruments and accessories

Hybrid instruments

HybridKnife, T-type, full-stream, \varnothing 2.3 mm, length 1.9 m

with connecting plug international (3-Pin)



No. 20150-060

= 5

2.5 kVp

30150-116



STERILE EO

HybridKnife, I-type, full-stream, \varnothing 2.3 mm, length 1.9 m

with connecting plug international (3-Pin)



No. 20150-061

= 5

2.5 kVp

30150-116



STERILE EO

HybridKnife, O-type, full-stream, \varnothing 2.3 mm, length 1.9 m

with connecting plug international (3-Pin)



No. 20150-062

= 5

0.99 kVp

30150-116



STERILE EO

Instruments and accessories

Hybrid instruments

HybridAPC, \varnothing 2.3 mm, length 1.9 m



No. 20150-015



= 5

4.3 kVp



30150-147



STERILE EO

Applicator, straight, with monopolar electrosurgical function, \varnothing 6 mm, length 80 mm
with suction, e.g. for liver surgery



No. 20150-036



= 4

2.5 kVp



30150-205



STERILE EO

Instruments and accessories

Applicators and probes


Applicator, curved tip, ø 5 mm, length 183 mm

e.g. for urology



No. 20150-025



 30150-000



STERILE EO


Applicator, curved tip, ø 5 mm, length 336 mm

e.g. for laparoscopy



No. 20150-026



 30150-000



STERILE EO


Applicator, straight, ø 6 mm, length 65 mm

with suction, e.g. for liver surgery



No. 20150-030



 30150-000



STERILE EO


Appikator, flexible tip, ø 6 mm, length 65 mm

with suction, e.g. for liver surgery



No. 20150-031



 30150-000



STERILE EO


Applicator, straight, ø 6 mm, length 306 mm

with suction, e.g. for TME



No. 20150-038



 30150-000



STERILE EO

Instruments and accessories

Applicators and probes

Applicator, straight, \varnothing 6 mm, length 180 mm

with suction, e.g. for TEM



No. 20150-039

= 5

30150-000



STERILE EO

Applicator, bayonet, \varnothing 2.8 mm, length 90 mm

with suction



No. 20150-041

= 5

30150-000



STERILE EO

Flexible probe, \varnothing 1.3 mm, length 2.2 m

e.g. for flexible endoscopy



No. 20150-020

= 5

30150-000



STERILE EO

Instruments and accessories

Accessories for ERBEJET 2

Pump (single-use)



No. 20150-300

= 5

30150-055

E0

Pump cartridge plus for ERBEJET 2



No. 20150-301

= 5

30150-400

STERILE E0

ERBEJET 2 one-pedal foot switch with ReMode, AP & IP X8 equipment



No. 20150-101

= 1

30150-100

max 75 °C



ERBEJET 2 two-pedal foot switch with ReMode, AP & IP X8 equipment



No. 20150-100

= 1

30150-100

max 75 °C



Further accessories for foot switches	Piece	No.
Hoop guard for VIO two-pedal foot switch	1	20189-306
Cable for ERBEJET 2 foot switch, length 5 m	1	20150-104

Power cord	Country	Length	Piece	No.
Power cord 10A	Germany	0.5 m	1	51704-046
Power cord 10A	Germany	2.5 m	1	51704-037
Power cord 10A	Germany	4 m	1	51704-032
Power cord 10A	Germany	5 m	1	51704-042
Power cord 10A	Great Britain	2.5 m	1	51704-038
Power cord 10A	Great Britain	5 m	1	51704-043
Power cord 10A	Italy	2.5 m	1	51704-040
Power cord 10A	Italy	5 m	1	51704-045
Power cord 10A	Canada	4 m	1	51704-039
Power cord 10A	India	4 m	1	20128-004

Instruments and accessories

Accessories for ESM 2

Suction container, 1.5 liters



for use with suction bag 20340-012

No. 20340-102

 = 1

Suction bag, 1.5 liters



for use with suction container 20340-102

No. 20340-012

 = 30



Suction bag, 2.5 liters



for use with suction container 20340-103

No. 20340-013

 = 30



Instruments and accessories

Accessories for ESM 2

ESM membrane filter



for use with silicone hose 20340-004

No. 20340-100

= 10



Silicone hose, length 600 mm



No. 20340-004

= 1

max 95 °C

max 134 °C

Rail bracket, plastic



for use with suction container 20340-102 and 20340-103

No. 20340-002

= 1

Switching valve, rail bracket and 2 silicone hoses, length 300 mm



No. 20340-003

= 1

max 95 °C

max 134 °C

Connecting cable

ECB connecting cable for ERBEJET 2, length 0.45 m

No.

51709-008

Fastening set

Fastening set EIP 2 / ESM 2 to VIO CART 20180-000

No.

20180-135

Instruments and accessories

Carts and accessories

VIO CART, system carrier for VIO family and modules



dimensions: 630 x 940 x 650 mm (WxHxD)
 weight: 31 kg
 4 conductive castors \varnothing 100 mm, with fixing brake
 1 power input
 3 power outlets

incl. following accessories:
 1 grab handle out of stainless steel
 2 cable holders
 1 foot switch bracket (for currently available foot switches)

No. 20180-000

= 1

80180-018



Hook for irrigation liquid bag on bracket / VIO Cart



No. 20325-002

= 1

max 95 °C

max 134 °C

Upper bracket for VIO CART, length 260 mm




No. 20180-040

= 1



Lower bracket for VIO CART, length 390 mm

for use with secretion suction



No. 20180-042

= 1



Fastening sets

Fastening set VIO D/S / APC 2 / ERBEJET 2 on VIO CART 20180-000

No. 20180-131

Fastening set EIP 2 / ESM 2 to VIO CART 20180-000

No. 20180-135

Fastening set on ERBEJET 2, connection to VIO / APC 2

No. 20150-050

Important information

The information, recommendations and other details contained in this medium ("information") reflect our state of knowledge and the state of the art in science and technology at the time the medium was created. The information is of a general nature, non-binding and intended solely for general information purposes; it does not constitute instructions for use.

The information and recommendations contained in this medium do not create any legal obligations on the part of Erbe Elektromedizin GmbH and its affiliated companies ("Erbe") or any other claims against Erbe. The information does not constitute any guarantees or other information on the nature of the product; such information requires an express contractual agreement with Erbe in each individual case.

Erbe is not liable for any kind of damage resulting from following the information on this medium, regardless of any legal grounds for liability.

Every user of an Erbe product is responsible for testing the respective Erbe product in advance with regard to its properties and its suitability for the intended type of use or purpose. Suitable type of use of the respective Erbe product can be found in the instructions for use for the respective Erbe product. The user is obliged to ascertain whether the instructions for use available to him/her correspond to the current status of the specific Erbe product. Devices may only be used in accordance with the instructions for use.

Information on settings, application sites, application duration and use of the respective Erbe product is based on the clinical experience of individual centers, physicians and clinical literature independent of Erbe as well as post market surveillance and empirical evidence. They constitute reference values, or Experienced Starting Settings, which must be checked by the user to ensure that they are suitable for the specifically planned application. Depending on the circumstances of a specific application, it may be necessary to deviate from the information provided. The user is responsible for checking this when using Erbe products. We wish to note that science and technology are subject to constant development as a result of research and clinical experience. Similarly, medicine and clinical application are constantly evolving as a result of research and technological development. This may also make it necessary for the user to deviate from the information provided on this medium.

This medium contains information about Erbe products that may not be approved/cleared in certain countries. The user of the respective Erbe product is obliged to inform himself/herself as to whether the Erbe product he/she is using is legally approved/cleared in his/her country and/or to what extent, if any, there are legal requirements for or restrictions on its use.

This medium is not intended for use in the USA.

Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Germany

Phone +49 7071 755-0
info@erbe-med.com
erbe-med.com
medical-videos.com